

比色槽 (Cell / Cuvette)

E

比色槽

圖片名稱 Picture and Name	外觀尺寸 (mm) HxWxD	光徑長 (mm) Path Length	光徑寬 (mm) Inside Width	容量 (ml) Volume	產品編號 (Product Number)		
					遠紫外光線 石英玻璃 (190-2500nm)	紅外光線 石英玻璃 (250-3500nm)	光學玻璃 (340-2500nm)
 白壁微量比色槽 (含PTFE塞)	45x12.5x7.5	5	5	0.85	Q-393	I-393	G-393
	45x12.5x12.5	10	5	1.7	Q-394	I-394	G-394
	45x12.5x22.5	20	5	3.5	Q-395	I-395	G-395
	45x12.5x32.5	30	5	5.2	Q-396	I-396	G-396
	45x12.5x42.5	40	5	7.0	Q-397	I-397	G-397
	45x12.5x52.5	50	5	8.75	Q-398	I-398	G-398
	45x12.5x102.5	100	5	17.5	Q-399	I-399	G-399

圖片名稱 Picture and Name	外觀尺寸 (mm) HxWxD	光徑長 (mm) Path Length	光徑寬 (mm) Inside Width	容量 (ml) Volume	產品編號 (Product Number)		
					遠紫外光線 石英玻璃 (190-2500nm)	紅外光線 石英玻璃 (250-3500nm)	光學玻璃 (340-2500nm)
 微量螢光比色槽 (含蓋)	45x12.5x12.5	10	1	0.35	Q-504	I-504	G-504
	45x12.5x12.5	10	2	0.70	Q-514	I-514	G-514
	45x12.5x12.5	10	3	1.00	Q-524	I-524	G-524
	45x12.5x12.5	10	4	1.40	Q-534	I-534	G-534
	45x12.5x12.5	10	5	1.70	Q-544	I-544	G-544

圖片名稱 Picture and Name	外觀尺寸 (mm) HxWxD	光徑長 (mm) Path Length	光徑寬 (mm) Inside Width	容量 (ml) Volume	產品編號 (Product Number)		
					遠紫外光線 石英玻璃 (190-2500nm)	紅外光線 石英玻璃 (250-3500nm)	光學玻璃 (340-2500nm)
 微量螢光比色槽 (含PTFE塞)	45x12.5x12.5	10	1	0.35	Q-554	I-554	G-554
	45x12.5x12.5	10	2	0.70	Q-564	I-564	G-564
	45x12.5x12.5	10	3	1.00	Q-574	I-574	G-574
	45x12.5x12.5	10	4	1.40	Q-584	I-584	G-584
	45x12.5x12.5	10	5	1.70	Q-594	I-594	G-594

圖片名稱 Picture and Name	外觀尺寸 (mm) HxWxD	光徑長 (mm) Path Length	光徑寬 (mm) Inside Width	容量 (ml) Volume	產品編號 (Product Number)		
					遠紫外光線 石英玻璃 (190-2500nm)	紅外光線 石英玻璃 (250-3500nm)	光學玻璃 (340-2500nm)
 壁光流動比色槽	45x12.5x7.5	5	Ø3	0.035	Q-63	I-63	G-63
	45x12.5x12.5	10	Ø3	0.07	Q-64	I-64	G-64
	45x12.5x22.5	20	Ø3	0.14	Q-65	I-65	G-65
	45x12.5x32.5	30	Ø3	0.21	Q-66	I-66	G-66
	45x12.5x7.5	5	Ø2	0.015	Q-67	I-67	G-67
	45x12.5x12.5	10	Ø2	0.03	Q-68	I-68	G-68
	45x12.5x22.5	20	Ø2	0.06	Q-69	I-69	G-69
	45x12.5x12.5	10	Ø1.5	0.018	Q-70	I-70	G-70