

圖片名稱 Picture and Name	外觀尺寸 (mm) HxWxD	光徑長 (mm) Path Length	光徑寬 (mm) Inside Width	容量 (ml) Volume	產品編號 (Product Number)		
					遠紫外光線 石英玻璃 (190-2500nm)	紅外光線 石英玻璃 (250-3500nm)	光學玻璃 (340-2500nm)
 黑壁半微量比色槽(含蓋)	45x12.5x7.5	5	4	0.7	Q-23	I-23	G-23
	45x12.5x12.5	10	4	1.4	Q-24	I-24	G-24
	45x12.5x22.5	20	4	2.8	Q-25	I-25	G-25
	45x12.5x32.5	30	4	4.2	Q-26	I-26	G-26
	45x12.5x42.5	40	4	5.6	Q-27	I-27	G-27
	45x12.5x52.5	50	4	7.0	Q-28	I-28	G-28
	45x12.5x102.5	100	4	14.0	Q-29	I-29	G-29

圖片名稱 Picture and Name	外觀尺寸 (mm) HxWxD	光徑長 (mm) Path Length	光徑寬 (mm) Inside Width	容量 (ml) Volume	產品編號 (Product Number)		
					遠紫外光線 石英玻璃 (190-2500nm)	紅外光線 石英玻璃 (250-3500nm)	光學玻璃 (340-2500nm)
 黑壁半微量比色槽(含PTFE塞)	45x12.5x7.5	5	4	0.7	Q-33	I-33	G-33
	45x12.5x12.5	10	4	1.4	Q-34	I-34	G-34
	45x12.5x22.5	20	4	2.8	Q-35	I-35	G-35
	45x12.5x32.5	30	4	4.2	Q-36	I-36	G-36
	45x12.5x42.5	40	4	5.6	Q-37	I-37	G-37
	45x12.5x52.5	50	4	7.0	Q-38	I-38	G-38
	45x12.5x102.5	100	4	14.0	Q-39	I-39	G-39

圖片名稱 Picture and Name	外觀尺寸 (mm) HxWxD	光徑長 (mm) Path Length	光徑寬 (mm) Inside Width	容量 (ml) Volume	產品編號 (Product Number)		
					遠紫外光線 石英玻璃 (190-2500nm)	紅外光線 石英玻璃 (250-3500nm)	光學玻璃 (340-2500nm)
 黑壁微量比色槽(含蓋)	45x12.5x7.5	5	2	0.35	Q-43	I-43	G-43
	45x12.5x12.5	10	2	0.7	Q-44	I-44	G-44
	45x12.5x22.5	20	2	1.4	Q-45	I-45	G-45
	45x12.5x32.5	30	2	2.1	Q-46	I-46	G-46
	45x12.5x42.5	40	2	2.8	Q-47	I-47	G-47
	45x12.5x52.5	50	2	3.5	Q-48	I-48	G-48
	45x12.5x102.5	100	2	7.0	Q-49	I-49	G-49

圖片名稱 Picture and Name	外觀尺寸 (mm) HxWxD	光徑長 (mm) Path Length	光徑寬 (mm) Inside Width	容量 (ml) Volume	產品編號 (Product Number)		
					遠紫外光線 石英玻璃 (190-2500nm)	紅外光線 石英玻璃 (250-3500nm)	光學玻璃 (340-2500nm)
 黑壁微量比色槽(含PTFE塞)	45x12.5x7.5	5	2	0.35	Q-53	I-53	G-53
	45x12.5x12.5	10	2	0.7	Q-54	I-54	G-54
	45x12.5x22.5	20	2	1.4	Q-55	I-55	G-55
	45x12.5x32.5	30	2	2.1	Q-56	I-56	G-56
	45x12.5x42.5	40	2	2.8	Q-57	I-57	G-57
	45x12.5x52.5	50	2	3.5	Q-58	I-58	G-58
	45x12.5x102.5	100	2	7.0	Q-59	I-59	G-59