



Eco 系列

一般型 (不具校正功能)



Airless 系列

適用有機溶劑



產品規格

LABMAX Eco 系列							
產品貨號	WTG-5370001	WTG-5370002	WTG-5370003	WTG-5370004	WTG-5370005	WTG-5370006	WTG-5370007OE
分配體積	0.25 ~ 2.5 mL	0.5 ~ 5.0 mL	1.0 ~ 10.0 mL	2.5 ~ 25.0 mL	5.0 ~ 50.0 mL	10.0 ~ 100.0 mL	25.0 ~ 250.0 mL
最小刻度	0.05 mL	0.1 mL	0.2 mL	0.5 mL	1 mL	2 mL	5 mL
刻度顯示	mL						
誤差值	系統誤差 (α)	± 0.012 mL	± 0.03 mL	± 0.06 mL	± 0.15 mL	± 0.3 mL	± 0.5 mL
	隨機誤差 (CV)	≤ 0.002 mL	≤ 0.005 mL	≤ 0.01 mL	≤ 0.025 mL	≤ 0.05 mL	≤ 0.1 mL
螺牙	GL32	GL32	GL32	GL45	GL45	GL45	GL45
隨附轉接頭	A25, A28, A45, S40			A25, A28, A32/45, S40			A25, A28, A32/45, A38, S40
高溫高壓滅菌	Yes	Yes	Yes	Yes	Yes	Yes	Yes

LABMAX Airless 系列						
產品貨號	WTG-5370801	WTG-5370802	WTG-5370803	WTG-5370804	WTG-5370805	WTG-5370806
分配體積	0.25 ~ 2.5 mL	0.5 ~ 5.0 mL	1.0 ~ 10.0 mL	2.5 ~ 25.0 mL	5.0 ~ 50.0 mL	10.0 ~ 100.0 mL
最小刻度	0.05 mL	0.1 mL	0.2 mL	0.5 mL	1 mL	2 mL
刻度顯示	mL					
誤差值	系統誤差 (α)	± 0.012 mL	± 0.03 mL	± 0.06 mL	± 0.15 mL	± 0.3 mL
	隨機誤差 (CV)	≤ 0.002 mL	≤ 0.005 mL	≤ 0.01 mL	≤ 0.025 mL	≤ 0.05 mL
螺牙	GL32	GL32	GL32	GL45	GL45	GL45
隨附轉接頭	A25, A28, A45, S40			A25, A28, A32/45, S40		
高溫高壓滅菌	Yes	Yes	Yes	Yes	Yes	Yes